

IMMUNOSTAINING PACKAGE

PRODUCT DESCRIPTION

P/N: O-SEQ-PCK

Version Jul. 24


This package includes all the components for sequential injections into a cell chamber or a chip.

Fluigent's product allow for an automated and precise perfusion with a short switching time. Up to 10 different samples can be put through the chip, followed by a sorting between waste and valuable product.



PLATFORM DESCRIPTION

- 1 LineUp Flow EZ pressure controller (1000 mbar)
- 1 LineUp LINK Module (software control)
- 1 FLOW UNIT M
- 1 2-SWITCH
- 1 M-SWITCH
- 1 SWITCH EZ
- 12 2-CAP series 15 mL
- 1 10 Position Pressure Manifold



TECHNICAL SPECIFICATIONS

Fluid handling system

Product	Part number
1*LINK module	LU-LNK-0002
1*Flow EZ 1000 mbar	LU-FEZ-1000
1*2-SWITCH	2SW003
1*SWITCH EZ	LU-LNK-0002
1*M-SWITCH	ESSMSW003

Reservoirs

Product	Part number
12*15 mL Pcap with 15 mL Falcon tube	P-CAP15-HP
1*10 Position Pressure Manifold	CTQ-MANI

Flow meters

Product	Part number
1*Flow unit M	FLU-M-D

Tubing

Product	Part number
12*Tubing and connection kit P-CAP 15 mL FEP tubing with an ID of 500 microns	CTQ-KIT-PC15

